

2040-ASTM, ADAPTER, 105° BEND PLAIN

Chemical Resistant

[Read More](#)

SKU:

Price:

Categories: [Adapters](#), [Laboratory Glassware](#)

Tags: [Adapter](#), [connector](#), [Lab Glassware](#), [Laboratory Glassware](#)



Product Description

Chemical Resistant

Part No.	FEMALE Joint	Pack QTY.
2040-14F	14/20	10
2040-19F	19/22	10
2040-24F	24/40	10
2040-29F	29/42	10

Features:

- Connecting Glassware:** The adapter may be used to connect two pieces of glassware or tubing at a specific angle.
- Adjusting Setup Orientation:** It could be used to alter the orientation of a particular part of a laboratory apparatus to optimize the experimental setup.
- Creating Custom Configurations:** Researchers often need to customize their setups for specific experiments, and adapters with specific bends can be valuable in achieving the desired configurations.
- Facilitating Fluid Flow:** In setups involving fluid transfer, a 105° bend adapter could be used to redirect the flow of liquids or gases.