



2050-ASTM, ADAPTER, 105° BEND

Chemical Resistant

[Read More](#)

SKU:

Price:

Categories: [Adapters](#), [Laboratory Glassware](#)

Tags: [Adapter](#), [connector](#), [Lab Glassware](#), [Laboratory Glassware](#)

Product Description

Chemical Resistant

Part No.	MALE & FEMALE Joint	Pack QTY.
2050-14FM	14/20	10
2050-19FM	19/22	10
2050-24FM	24/40	10
2050-29FM	29/42	10

Applications for an adapter with a 105° bend in a laboratory setting:

- Customized Experimental Setups:** Adapters with specific bends are often used in laboratories to create customized experimental setups. Researchers may need to achieve specific angles or orientations in their apparatus, and adapters like this could facilitate such configurations.
- Glassware Connection:** Adapters are commonly used to connect different pieces of laboratory glassware. A 105° bend might be useful in creating a connection between two components at a non-standard angle.
- Fluid Transfer:** In setups involving the transfer of liquids or gases, adapters with bends can be used to redirect the flow in a controlled manner. The 105° bend could be part of a fluid-handling system.
- Space Optimization:** Sometimes, laboratory setups need to be compact or require a specific spatial arrangement. Adapters with bends can help optimize space and allow for more efficient use of the laboratory workspace.